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## **Instruction for use of DURASPLINT<sup>®</sup> LC**

The transparent light-curing resin  
specifically designed to adjust  
DURAN<sup>®</sup> splints.

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We recommend the above mentioned items for the fabrication of adjusted occlusal splints made of DURAN® to facilitate your work and to ensure a user-friendly production process.

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### Instructions regarding work safety

Avoid skin contact with uncured resins. Make sure to wear gloves during processing.


Wear protective clothing (smock, protective goggles, safety mask). Switch on suction unit.

### Instructions regarding light polymerisation

Light polymerisation requires a wavelength between 350 nm and 450 nm. The LC-6 Light Oven offers this spectral range.

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Use flat trimmed models. Mount models in an articulator following the dentist's instructions and using the recommended bite registration. Optionnally, you can use magnetic systems.

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Measure the model using the parallelometer taking into account the direction of insertion and mark the prosthetic equator.

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Mark splint extension in the area of slight undercuts in red.

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We recommend blocking out strong undercuts, like space closures, using BLUE-BLOKKER® or SIL-KITT to facilitate removing the finished splint.

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The same applies for interdental spaces and bridge elements.

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Pressurise an ISOFOILAN® foil first before thermoforming DURAN® 1.5 mm or 2.0 mm according to the instructions.

We recommend using the pellet cover.

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**Important:** It is essential to use DURAN® as basis for the light-curing material.

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Reduce the pressure-moulded DURAN® roughly by means of the HM carbide bur or the diamond separating disc.

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**Hint:** For optimum fit remove splint only after light-polymerisation of DURASPLINT® LC and after reducing the splint.

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The material surface in the complete splint extension should be sandblasted using aluminium oxide ( $Al_2O_3$ ) or roughened with a DIMO® wheel (brown, coarse).

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Take off grease residues from roughened area with ethanol to guarantee a good compound. Apply one thin layer of LC-Primer on the complete splint extension and proceed with curing for 5 minutes in LC-6 Light Oven.

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Apply LC Insulation to the antagonist model with a brush, to prevent DURASPLINT® LC material from sticking to the model.

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Adapt one bar of the DURASPLINT® material on the pressure moulded splint.

Additionally, CLEAR-BLOKKER® can be applied as modelling aid to reach a smooth result.

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**Important:** Avoid skin contact with uncured material. Wear gloves while processing!

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For easier modelling moisten your fingertips with LC Insulation. Now you can model the material to the desired shape.

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**Hint:** We recommend applying a thin layer of the material, slightly exceeding the splint extension, in order to avoid visible material transitions between DURASPLINT® LC and DURAN®.

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Insert movement patterns and impressions using the articulator. Movements should be made only once in the corresponding direction before returning to the centric position with the articulator being open.

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**Hint:** Repeated movements or return movements with the articulator being closed might cause distortion due to the slightly sticky material.

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Remove the model with the pressure moulded splint from articulator and place it centrally in the LC-6 Light Oven. Then cure the material 2 times for 10 minutes. If the LC-6 Light Oven is used for a serial workflow, please make sure to cool down the LC-6 Light Oven after each polymerisation period by opening the drawer.

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After light-curing you can finish the splint still being on the working model.

The slightly sticky inhibition layer can be removed by using medical ethanol (Isopropyl).

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
For initial finishing we recommend using LC Burs A, B or C

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After grinding in sliding surfaces and cusps according to the therapeutic instructions you can proceed with finishing the splint.

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Reduce material up to the splint extension carefully by means of the HM Carbide Cutter or the diamond separating disc. Avoid any contact between carbide cutter and model.

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Remove splint carefully from the model, taking it off simultaneously from both sides.

**Hint:** This is essential to avoid cracks that might occur if the splint is removed from one side only.

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A chamfered POLYFIX disc is best suited for interdental finishing of the splint.

For further finishing we recommend using the LC burs A, B and C.

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We recommend using the DIMO® disc coarse to smooth the surfaces of the rims or optionally sandpaper (both items are not included in the DURASPLINT® LC kit).

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Polished DURAN® splint with DURASPLINT® LC build-up on DURAN®.

**Hint:** We do not recommend occlusal polishing, in order to maintain the full functionality of the adjusted splint surfaces.

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You can use CLEAR-BLOKKER® in a layer of 1-2 mm for any kind of completions or corrections. Make sure the inhibition layer stays on the object. If this layer has already been removed, apply LC-Primer to the corresponding area and cure for 5 minutes. Then apply CLEAR-BLOKKER® directly from the syringe using an application cannula. Proceed with light-curing for 10 minutes, e. g. in the LC-6 Light Oven.

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### Additional information:

DURASPLINT® LC and CLEAR-BLOKKER® do not contain methyl methacrylate monomer (only the LC-Primer contains methyl methacrylate monomer). Please use only light curing ovens that do not develop high temperature levels during operation (we recommend the LC-6 Light Oven).

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The CETRON® Care Set is the best choice for daily cleansing and care.

CE