

# Finishing instruments

## BIOPERM-Trimmer / POLYFIX

Trimming and initial polishing of soft-elastic and sandwich material.

## OSAMU polisher

Finishing and polishing of soft-elastic and sandwich material.

### Delivery program

<b>BIOPERM trimmer, blue</b>	2 pieces	<b>#3226</b>
<b>POLYFIX on mandrel, brown</b>	1 piece	<b>#3371</b>
<b>POLYFIX without mandrel, brown</b>	6 pieces	<b>#3372</b>
<b>OSAMU polisher, grey</b>	2 pieces	<b>#3247</b>



## DIMO®

Soft, flexible abrasive wheels with open non woven fibre structure for professional finishing of pressure moulding material, hard and soft resins and metals. They are equally perfectly suited for grinding and smoothing of silicone, polyurethane and permanently soft materials. Optimized material, minimized unravelling and low thermal development guarantee a much longer life of the wheels.

### Application

Finishing and initial polishing of elastomers, sandwich materials, mouth guard and orthodontic appliances.

### Delivery program

<b>DIMO® brown, coarse</b>	10/50 pieces	<b>#3382</b>
<b>DIMO® green, medium</b>	10/50 pieces	<b>#3383</b>
<b>DIMO® maroon, fine</b>	10/50 pieces	<b>#3384</b>
<b>DIMO® mandrel</b>	10 pieces	<b>#3380</b>

**The all-rounders  
with longer  
endurance**



## DIMO®PRO / DIMO®PRO SLIM

It's the combination of a special elastomer compound and fine fibre structure that makes DIMO®PRO so versatile. You can achieve the best results in grinding and smoothing as well as in polishing using different speeds.

### Application

**DIMO®PRO:** Finishing of large-scale marginal areas and for final and high-gloss polishing of splints, positioners, mouth guard and orthodontic appliances. Provides highly-polished rims and uniform surfaces of hard/soft sandwich material.


**DIMO®PRO SLIM:** With tapered rim for finishing splints, best suited for the gaps in the area of the frenulum and for interdental finishing.

### Delivery program

<b>DIMO®PRO</b> Wheel 24 × 6 mm	10/50 pieces	<b>#3381</b>
<b>DIMO®PRO SLIM</b> Wheel 24 × 6 mm	10/50 pieces	<b>#3376</b>

## Finishing set

Set comes with eight grinding wheels for cutting and finishing of soft elastic and hard elastic material. With color coding for easy identification.



### // Finishing set

1 piece

#3378

Set comes with one each cutting bur, finishing bur, HM carbide bur, HM tricutter fine, BIOPERM trimmer, POLYFIX, OSAMU polisher, DIMO®PRO SLIM and dressing stone.

## Cutting bur

Spiral cutting bur for quick, neat cutting of hard pressure moulding material.

## HM Carbide cutter

Triangular carbide cutting bur for quick, neat cutting of hard pressure moulding material.

## HM Tricutter fine

Fine triangular cutter for interdental finishing of splints.

## Finishing bur

Square cutting bur made of special soft alloy that will not damage clasps and wires.

## Delivery program

// Cutting bur	6 pieces	#3214
// HM Carbide cutter, yellow	6 pieces	#3369
// HM Tricutter fine, green	6 pieces	#3370
// Finishing bur, red	6 pieces	#3377

