

STEP 1 – BULK COMPOSITE REDUCTION



- 18 fluted carbide
- Quick removal of bulk composite
- Gentle on enamel
- Use with or without H2O



STEP 2 – RESIDUAL COMPOSITE REMOVAL



- Polishes tooth surface
- Black point discloses when composite is fully removed
- Removes residual composite
- Does not create heat



STEP 3 – FINISHING & POLISHING



- Pre-packaged—sterile cups for individual patient use. Pre-mixed—just open and apply with prophyl cup.
- No added chemicals, flavoring, fluoride or oils.
- Enamel polisher—1st & Final[®] can be used as final polishing step with the Reliance Renew[™] Finishing System.
- Restore Diamond Polishing Paste for restoration of glaze on artificial surfaces.

